890 Series – 3A Sanitary Pressure Transducer

- ▶ For Clean-In-Place (CIP) and Sterilize-In-Place (SIP)
- ▶ 0.20% Full Scale Accuracy
- ▶ No Liquid Fill Diaphragms

The 890 Series meets 3A sanitary design standards and is fully sealed to withstand external high pressure washdowns. These units are packaged in rugged welded stainless steel housings and are exceptionally insensitive to vibration, shock and environmental extremes. A small size and tri-clover sanitary pressure fitting allow direct mounting in most CIP and SIP installations. Other features include IC-based circuitry, a 1/2" NPT conduit fitting and shielded cable with vent tube. Sealed screws provide access to zero and span adjustments.

Specifications

•	
Input	
Pressure Range	1 to 1000 psig
Proof Pressure	see ordering chart
Burst Pressure	see ordering chart
Fatigue Life	>1 million cycles
Performance	
Output	4-20 mA (2 Wire)
Supply Voltage (Vs)	18-38 Vdc
Accuracy	0.20% FS
Thermal Error Zero	.02% FS/°F (.036%FS/°C)
Thermal Error Span	.02% FS/°F (.036%FS/°C)
Compensated Temperatures	20°F to 180°F (-7°C to +80°C)
Operating Temperatures	-40°F to +260°F (-40°C to +125°C)
Storage Temperatures	-65°F to +260°F (-54°C to +127°C)
Zero Tolerance	1% FS (±0.5 mA adjustable)
Span Tolerance	1% FS (±0.5 mA adjustable)
Maximum Loop Resistance	(Vs-18) x 50
Response Time	10 ms
Mounting Effects	.15% FS (.25% FS for 1.5" Tri-Clover)
Mechanical Configuration	
Pressure Port	1.5" or 2" Tri-Clover Sanitary Fitting
Wetted Parts	316 Stainless Steel
Electrical Connection	1/2" NPT Conduit Fitting and Strain Relief with 15 ft. Cable
Enclosure	Stainless Steel
Vibration	10g Peak Sinusoidal, 50 to 1000 Hz
Acceleration	10g
Shock	50g
Approvals	Meets 3-A Sanitary Standards
Weight	8 oz



Applications

- · Food Processing
- · Dairy & Beverage Processing
- Pharmaceutical Processing
- Sanitary Pipelines


How They Operate

A stainless steel diaphragm and an insulated electrode form a variable capacitor. Pressure on the diaphragm alters the sensor's capacitance, which is then detected and converted to a highly accurate linear 4-20 mA signal by electronic circuitry featuring Gems' patented charge-balance principle. Low hysteresis, very stable operation and negligible clamping effect are inherent.



Dimensions





2"Fitting

Gems adheres to strict quality standards including MIL-1-45208A and ANSI-2540-1.

How to Order

Order as 890 Series Sanitary Pressure Transmitters. Specify Pressure Range (tablulated below), Fitting Size and any Options. Use bold characters to construct a product code.

C890 - 10 - 1.5 - 884 -822

SELECT

1. Series — **C890** - 890 Series

2. Pressure Ranges

1.5" Tri-Clover Sanitary Fittings 2" Tri-Clover Sanitary Fittings **Operating Range** Proof Burst **Operating Range** Proof **Burst** psig in. H₂0 psig psig psig psig psig 1 22.7 50 100 30 1000 1200 2 55.4 100 150 60 1000 1200 5 200 100 1000 138.4 150 1200 10 276.8 150 200 300 1000 1200 15 415.2 150 200 500 1000 1500 30 1000 2400 830.4 150 300 1250 60 1160.8 180 400 -14.7 to 15 1000 1200 100 -14.7 to 45 2768.0 200 400 1000 1200 150 4152.0 225 400

300

3. Pressure Port

1.5 - 1.5" Tri-Clover Sanitary Fitting

150

2.0 - 2" Tri-Clover Sanitary Fitting

-14.7 to 15 -407 to 415

4. Optional

715 - ±0-1% FS (RSS) accuracy

884 - 20 Ra finish

911 - Etched metal stainless steel tag

5. Cable Length

816-825 - For cable lengths of 16 to 25 feet (15 ft. is standard).
Please specify cable length by code (e.g., 820 for 20 ft. cable).
Consult factory for cable longer than 25 feet.